



Use Able AI chat for product recommendations

Chat now



Q Enter a keyword, multiple keywords, or ask a question.



The 125 mL and 500 mL sizes are provided in narrow-neck bottles.

[View All Images \(11\)](#)



H&E Staining Kit (Hematoxylin and Eosin)

Cat No. PK10031

Publications

8

[Review Product](#)

Subject

Tissue sections and tissue cells

Specificity

Non-toxic, produces clear staining results, simple, quick and ready to use

Size: 10mL

10mL

125mL

500mL

\$55/_{10mL}

In stock in Manchester, UK. Delivery in 1 week.

Freight/Packing: \$40.00

Quantity

-

1

+

[Request bulk or custom quote](#)

[Add to Basket](#)

Proteintech Guarantee

Overview

[Product Information](#)[Storage](#)[Usage](#)[Cautions](#)[Documentation](#)[Publications \(8\)](#)[Print datasheet](#)

Product Information

Proteintech's PK10031 (H&E staining kit) can be used for tissue and cell staining. Cell nuclei can be stained blue while cytoplasm can be stained red. The formulas of the reagents in the kit have been optimized to avoid using toxic chemicals, e.g. mercury or methanol. Staining with the H&E Staining Kit is fast and easy.

The product is provided in a ready-to-use format. No further dilution is required. Each milliliter of this product can stain ~10 slides.

Component	Size		
Mayer's Hematoxylin Solution	10 mL	125 mL	500 mL
Eosin Stain Solution	10 mL	125 mL	500 mL
Bluing Buffer	10 mL	125 mL	500 mL

Storage

Ambient. Stable for one year.

Usage

For FFPE samples:

1. Immerse slides in xylene for 10 minutes. Repeat in fresh xylene again for another 10 minutes. Then rehydrate the slides by sequentially incubating in 100%, 95%, 80%, and 60% ethanol solutions for 5 minutes each. Finally wash the slides with distilled water for 2 minutes.
2. Apply Mayer's Hematoxylin Solution from the kit onto the tissue slides and make sure the tissues are fully covered. After staining for 1-2 minutes, wash the slides with distilled water. (If staining is too dark, the staining could be lightened by immersing the slides in 1% HCl ethanol solution). Apply

bluing reagent from the kit to the slides for 2 minutes. Once the sample turns blue, wash the slides briefly with distilled water 3 times.

3. Apply Eosin staining solution from the kit to the slides and make sure the tissues are fully covered. Stain for 2-5 minutes. Wash the slides with distilled water after staining.

4. Dehydrate the slides with gradient ethanol followed by clearing, mounting and observation.

For frozen samples:

1. Wash the slides with distilled water, apply Mayer's Hematoxylin Solution from the kit onto the tissue slides and make sure the tissues are fully covered. After staining for 1-2 minutes, wash the slides with distilled water. (If staining is too dark, the staining could be lightened by immersing the slides in 1% HCl ethanol solution). Apply bluing reagent from the kit to the slides. Blue for 2 minutes. Once the sample turns blue, wash the slides briefly with distilled water for 3 times.

2. Apply Eosin staining solution from the kit to the slides and make sure the tissues are fully covered. Stain for 2-5 minutes. Wash the slides with distilled water after staining.

3. Dehydrate the slides with gradient ethanol followed by clearing, mounting and observation.

Cautions

1. Stain with 1-2 samples to make sure the staining achieves optimal results.

2. Wear a lab coat and gloves during operation.

Cited in Article as

PK10031, H&E Staining Kit (Hematoxylin and Eosin), Proteintech, IL, USA



Documentation

SDS

[H&E Staining Kit \(Hematoxylin and Eosin\)](#)

DATASHEET

[H&E Staining Kit \(Hematoxylin and Eosin\) Datasheet](#)

Publications

APPLICATION	TITLE
	Front Bioeng Biotechnol Photocrosslinked dual-network hydrogel for sutureless corneal stromal lenticule Implantation. Authors - Xianglong Yi View Article
	Mol Neurobiol HINT2-Mediated Mitochondrial Modulation Contributes to G-CSF Neuroprotection in Alcohol Use Disorder-Related Ischemic Stroke. Authors - Zhi-Chen Guo

APPLICATION	TITLE
	View Article
J Bioenerg Biomembr	KunMingShanHaiTang formula reprograms macrophage metabolism and promotes M2 polarization via the HIF-1α pathway to alleviate ulcerative colitis symptoms in a rat model. Authors - Zhiyun Zhang View Article
J Pain Res	Effect of Pulsed Radiofrequency with Various Time-Course on Rheumatoid Arthritis Joint Pain in Rats. Authors - Xin Li View Article
Cell Div	Linalool targets NR3C2 to inhibit NF-κB-mediated gastric cancer progression. Authors - Ying Hu View Article
Mol Pharm	Dihydroquercetin Nanomicelles Mitigate Hippocampal Apoptosis and Alleviate Autism-Like Behaviors in ASD Rats via the IKK/IκB/NF-κB Signaling Pathway. Authors - Shan Zhao View Article

[See more citations](#)



[Contact Us](#) [Privacy Policy](#) [Terms of use](#) [Terms and Conditions](#) [Trademark Information](#) [Imprint \(Impressum\)](#) [Modern Slavery Statement](#) [ChromoTek GmbH Privacy Policy](#)

Copyright © 2002-2026 Proteintech Group, Inc. All rights reserved.